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STS - Subdermal Therapy System

Clinical trial

Medical Research Department
SOR Internacional, S.A.

SORISA® 

Introduction

Sorisa STS is a platform for the treatment of cellulite and body volume excess, combining several non-invasive, painless techniques that act at different levels, achieving significant, long-lasting results in silhouette remodelling.

The STS comprises:

- a) Vacuum system to drain liquid retention (Alster and Tanzi, 2005)
- b) 3 MHz Ultrasounds for speeding up tissue drainage (Siems et al. 2005)
- c) Red LEDs for cell activation and drainage (Sasaki et al. 2007)
- d) MMF (modulated medium frequency) currents to activate circulation and enhance product penetration by means of poration (Amat A, Corral-Baqués MI, 2007).

SORISA® 

Material and methods

The study of the STS equipment was performed at 5 collaborating centres at Barcelona, Tarragona and Madrid.

Sample size: 48 patients

Body measurements: To assess the results of the treatments in terms of circumference loss, the silhouettes of the areas treated were measured before and after treatment

Digital photographs: Digital pictures were obtained before and after the treatment for the clinical assessment of the result


Clinical study

Dermal scanner by 22 MHz ultrasounds (DUB System, from Taberna Pro Medicum, Germany): It analyses the degree of dermal compactation. Compact skin is seen as uniformly green. Skin with dark areas on the scanner image is not very compact and has a lower density. This may be due to a process such as cellulite, which causes accumulations of liquid and fat in the skin.

[H.C. Corting et al. 1999; P. Altmeyer et al. 1991; B. Mole et al. 2004]

Treatment Gel:

In all cases the G5 treatment gel from Sorisa was applied, which contains as most significant components: Laminaria digitata, paullinia cupana, carnitina, ruscus aureolatus, fucus vesiculosus, equisetum arvense, chondrus crispus, arnica montana, potassium iodide, hederia helix, panax ginseng, centella asiatica.

Clinical study


Treatment parameters:

- 1) Programme: the most appropriate cellulite programme for each subject was chosen (compact cellulite, flaccid, edematous, or mixed).
- 2) MFM currents:
 - 20 minutes per area
 - Current percentage: 23-65 %
- 3) Ultrasounds:
 - 7 minutes per area
 - 1-2 W/cm²
 - 3 MHz
- 4) Vacuum:
 - 14 minutes per area
 - Pressure: 145-365 mbar
- 5) Sessions: 2 or 3 sessions/week, depending on the severity of cellulite

Clinical study

Results

Reduction in centimeters of the abdominal perimeter, waist and around the femoral trochanter, after an average of 12 sessions treatment with STS equipment

Results

Photographic assessment

Treatment file
Patient: A 40-year old man
Sessions: 12
Programme: Body remodelling
Mean level of MFM currents: 60%
Vacuum pressure: 180-320 mbar (continuous, pulsed and undulated mode)
Application: Plates in abdomen
Abdominal perimeter -4 cm

Treatment file
Patient: A 46-year old woman
Sessions: 12
Programme: Compact cellulite
Mean level of MFM currents: 42%
Vacuum pressure: 145-260 mbar (continuous, pulsed and undulated mode)
Application: Plates in abdomen, buttocks and thigh
Perimeter of abdomen -11 cm Perimeter of waist -13 cm

Results Photographic assessment

SORISA® 

Treatment file

Patient: A 46-year old woman

Sessions: 12

Programme: Compact cellulite

Mean level of MFM currents: 42%

Vacuum pressure: 145-260 mbar (continuous, pulsed and undulated mode)

Application: Plates in abdomen, buttocks, and thighs.

Perimeter of abdomen -11 cm Perimeter of waist -13 cm

Treatment file

Patient: A 48-year old woman

Sessions: 12

Programme: Mixed cellulite

Mean level of MFM currents: 43%

Vacuum pressure: 180-300 mbar (continuous, pulsed and undulated mode)

Application: Plates in abdomen and buttocks.

Perimeter of abdomen -7 cm Perimeter of waist -7 cm

Clinical study

Treatment file

Patient: A 31-year old woman

Sessions: 12

Programme: Compact cellulite

Mean level of MFM currents: 28%

Vacuum pressure: 150-260 mbar (continuous, pulsed and undulated mode)

Application: Plates in abdomen, legs and buttocks.

Perimeter of abdomen -4 cm Perimeter of waist -7 cm

Treatment file

Patient: A 36-year old woman

Sessions: 12

Programme: Mixed cellulite

Mean level of MFM currents: 63%

Vacuum pressure: 140-290 mbar (continuous, pulsed and undulated mode)

Application: Plates in abdomen, buttocks, and legs.

Perimeter of abdomen -4 cm Perimeter of waist -5 cm

Results Photographic assessment

SORISA® 

Treatment file

Patient: A 36-year old woman

Sessions: 12

Programme: Mixed cellulite

Mean level of MFM currents: 63%

Vacuum pressure: 140-290 mbar (continuous, pulsed and undulated mode)

Application: Plates in abdomen, buttocks, and legs.

Perimeter of abdomen -4 cm Perimeter of waist -5 cm

Treatment file

Patient: A 36-year old woman

Sessions: 12

Programme: Mixed cellulite

Mean level of MFM currents: 63%

Vacuum pressure: 140-290 mbar (continuous, pulsed and undulated mode)

Application: Plates in abdomen, buttocks, and legs.

Perimeter of abdomen -4 cm Perimeter of waist -5 cm


Clinical study

SORISA® 


Results
Dermal Scanner 22 MHz


Treatment file
Patient: A 27-year old woman
Sessions: 12
Area: Buttocks


Treatment file
Patient: A 46-year old woman
Sessions: 12
Area: Buttocks

Clinical study

Results Dermal Scanner 22 MHz


Treatment file
Patient: A 46-year old woman
Sessions: 12
Area: Thigh


Treatment file
Patient: A 33-year old woman
Sessions: 12
Area: Abdominal

Clinical study


Treatment file
Patient: A 56-year old woman
Sessions: 12
Area: Buttocks


Treatment file
Patient: A 60-year old woman
Sessions: 12
Area: Abdominal

Clinical study

Results Dermal Scanner 22 MHz

Treatment file
Patient: A 36-year old woman
Sessions: 12
Area: Thighs

Treatment file
Patient: A 36-year old woman
Sessions: 12
Area: Buttocks

Clinical study

SORISA® 

Conclusions

- The STS equipment is a combined system of techniques that have been shown to be effective reducing localised cellulite and fat deposit, thus contributing significantly to remodelling body silhouette.
- It is a non-invasive, painless, safe technique that provides results that can be proven over a significantly short time.

Clinical study

SORISA® 

References

- Alster TS, Tanzi EL. Cellulite treatment using a novel combination radiofrequency, infrared light, and mechanical tissue manipulation device. *J Cosmet Laser Ther.* 2005 Jun;7(2):81-5.
- Siems W et al., Anti-fibrosclerotic effects of shock wave therapy in lipedema and cellulite, *Biofactors.* 2005;24(1-4):275-82
- Sasaki GH et al., The effectiveness and safety of topical PhotoActif phosphatidylcholine-based anti-cellulite gel and LED (red and near-infrared) light on Grade II-III thigh cellulite: a randomized, double-blinded study, *J Cosmet Laser Ther.* 2007 Jun;9(2):87-96
- Amat A, Corral Baqués MI, Estudio Clínico Multicéntrico PSP System, Sorisa-SPAIN, 2007)
- Korting H C et al., *Ultraschall in der Dermatologie, Ein Atlas*, 1999, Blackwell Wissenschafts-Verlag Berlin-Wien
- Altmeyer P et al., *Ultrasound in Dermatology*, 1991, Springer-Verlag, Berlin
- Mole B et al., High frequency ultrasonography and celluscore: an improvement in the objective evaluation of cellulite phenomenon, 2004, *Ann Chir Plas Esthet.* Aug; 49(4): 387-95

The logo for SORISA, featuring the word "SORISA" in a serif font followed by a stylized graphic of three vertical bars in green and blue.

Sor Internacional, S.A.
Moianès, 13 - Pol. Ind. Can Casablanas
08192 Sant Quirze del Vallès - Barcelona - España
Tel. 93 721 69 20 - Fax 93 721 47 74
e-mail: sorisa@sorisa.com
www.sorisa.com

This document is the property of Sor Internacional, S.A. Its disclosure or reproduction without express consent from its owner is prohibited.

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